



# HeartStart Home to Hospital

A complete portfolio of cardiac resuscitation products and services for virtually any caregiver in virtually any environment.



13.5 lbs.

## HeartStart MRx Monitor/Defibrillator

- \_\_\_\_\_ Motion-tolerant NBP and SpO<sub>2</sub>
- \_\_\_\_\_ Microstream® etCO<sub>2</sub>
- \_\_\_\_\_ IBP (2 channels)
- \_\_\_\_\_ Temperature monitoring
- \_\_\_\_\_ CPR measurement & feedback
- \_\_\_\_\_ 12-lead ECG – all waves viewable on-screen
- \_\_\_\_\_ ST/AR arrhythmia detection
- \_\_\_\_\_ Noninvasive pacing and cardioversion
- \_\_\_\_\_ Monitor network connectivity
- \_\_\_\_\_ <5-second charge time
- \_\_\_\_\_ Long-life battery 8 hours – 1 battery
- \_\_\_\_\_ 5-year warranty
- \_\_\_\_\_ Quick Shock



14.7 lbs.

## HeartStart XL+ Monitor/Defibrillator

- \_\_\_\_\_ Manual & AED mode
- \_\_\_\_\_ 3/5 lead monitoring
- \_\_\_\_\_ SpO<sub>2</sub> Monitoring
- \_\_\_\_\_ NIBP monitoring
- \_\_\_\_\_ Non-invasive pacing & cardioversion
- \_\_\_\_\_ Quick Shock
- \_\_\_\_\_ 3-hour minimum battery life
- \_\_\_\_\_ 5-year warranty



3.4 lbs.

## HeartStart FR3

- \_\_\_\_\_ Adult and pediatric defibrillation with infant/child key
- \_\_\_\_\_ Optional manual override
- \_\_\_\_\_ Extensive daily self-tests
- \_\_\_\_\_ 12 hours operating time – 1 battery
- \_\_\_\_\_ Pre-connected pads with pads testing
- \_\_\_\_\_ Wireless data transfer
- \_\_\_\_\_ Exceptional ruggedness
- \_\_\_\_\_ SMART CPR
- \_\_\_\_\_ Quick Shock
- \_\_\_\_\_ 5-year warranty



3.5 lbs.

## HeartStart FRx

- \_\_\_\_\_ Adult and pediatric defibrillation with infant/child key
- \_\_\_\_\_ Pre-connected pads with pads testing
- \_\_\_\_\_ Easy data transfer
- \_\_\_\_\_ Exceptional ruggedness
- \_\_\_\_\_ Extensive daily self-tests
- \_\_\_\_\_ CPR Coaching
- \_\_\_\_\_ 4-year battery (standby life)
- \_\_\_\_\_ 8-year warranty



3.3 lbs.

## HeartStart OnSite / HeartStart Home

- \_\_\_\_\_ Adult and pediatric defibrillation
- \_\_\_\_\_ Integrated pads cartridge with pads testing
- \_\_\_\_\_ Easy data transfer
- \_\_\_\_\_ Extensive daily self-tests
- \_\_\_\_\_ CPR coaching
- \_\_\_\_\_ Quick Shock
- \_\_\_\_\_ 4-year battery (standby life)
- \_\_\_\_\_ 8-year warranty

For more information about our resuscitation solutions, contact Philips at 1-800-453-6860.

[www.philips.com/heartstart](http://www.philips.com/heartstart)

All HeartStart defibrillators are AHA Guidelines 2010 compliant

**PHILIPS**  
sense and simplicity